

Safety Data Sheet

according to UK REACH (SI 2020/1577) as amended

Date of issue: 11.09.2026

Version number: 1.1 (replaces version 1.0)

Revision date: 07.11.2012

* SECTION 1: Identification of the substance/mixture and of the company/undertaking

- **1.1 Product identifier**
- **Trade name:** Fmoc-Gln(trt)-OH >98.5 %
- **Product Code:** 60-1132-25
- **CAS Number:**
132327-80-1
- **EC number:**
603-564-1
- **1.2 Relevant identified uses of the substance or mixture and uses advised against**
- **Product category** PC21 Laboratory chemicals
- **Application of the substance / the mixture** Laboratory chemicals
- **Uses advised against** None identified
- **1.3 Details of the supplier of the safety data sheet**
- **Supplier:**
Severn Biotech Ltd.
Unit 2,
Park Lane,
Kidderminster,
Worcestershire.
DY11 6TJ
UK
Tel: 0044 1562 825286
Fax: 0044 1562 825284
email: info@severnbiotech.com
- **Further information obtainable from:** Product safety department.
- **1.4 Emergency telephone number:**
Members of the public seeking specific information on poisons should contact:
In England and Wales: NHS 111 - dial 111
In Scotland: NHS 24 - dial 111

* SECTION 2: Hazards identification

- **2.1 Classification of the substance or mixture**
- **Classification according to GB-CLP**
Skin Sens. 1 H317 May cause an allergic skin reaction.
Aquatic Chronic 4 H413 May cause long lasting harmful effects to aquatic life.
- **Additional information:** The full text of the classification and labelling elements mentioned is given in Section 16.
- **2.2 Label elements**
- **Labelling according to GB-CLP** The substance is classified and labelled according to the GB CLP regulation.
- **Hazard pictograms**



GHS07

- **Signal word** Warning
- **Hazard statements**
H317 May cause an allergic skin reaction.
H413 May cause long lasting harmful effects to aquatic life.
- **Precautionary statements**
P261 Avoid breathing dust.
P273 Avoid release to the environment.
P280 Wear protective gloves/protective clothing/eye protection/face protection.

(Contd. on page 2)

GB

Safety Data Sheet

according to UK REACH (SI 2020/1577) as amended

Date of issue: 11.09.2026

Version number: 1.1 (replaces version 1.0)

Revision date: 07.11.2012

Trade name: Fmoc-Gln(trt)-OH >98.5%

(Contd. of page 1)

P302+P352 IF ON SKIN: Wash with plenty of soap and water.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

P501 Dispose of contents/container in accordance with local regulations.

- **2.3 Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable
- **vPvB:** Not applicable

SECTION 3: Composition/information on ingredients

- **3.1 Substances**
- **CAS No. Description**
CAS: 132327-80-1 N2-[(9H-Fluoren-9-ylmethoxy)carbonyl]-N-(triphenylmethyl)-L-glutamine
- **Identification number(s)**
- **EC number:** 603-564-1

* SECTION 4: First aid measures

- **4.1 Description of first aid measures**
- **General information:** Immediately remove any clothing soiled by the product.
- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:**
Immediately wash with water and soap and rinse thoroughly.
If skin irritation continues, consult a doctor.
- **After eye contact:**
Check for and remove any contact lenses.
Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor.
- **After swallowing:**
Rinse out mouth and then drink plenty of water.
If symptoms persist consult doctor.
- **Information for doctor:** Treat symptomatically and supportively.
- **4.2 Most important symptoms and effects, both acute and delayed** Allergic reactions
- **4.3 Indication of any immediate medical attention and special treatment needed**
No further relevant information available.

* SECTION 5: Firefighting measures

- **5.1 Extinguishing media**
- **Suitable extinguishing agents:**
CO₂, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
Use fire extinguishing methods suitable to surrounding conditions.
- **For safety reasons unsuitable extinguishing agents:** Water with full jet
- **5.2 Special hazards arising from the substance or mixture**
As with most organic solids, fire is possible at elevated temperatures or by contact with an ignition source.
Fine dust dispersed in air in sufficient concentrations, and in the presence of an ignition source is a potential dust explosion hazard.
In case of fire, the following can be released:
Carbon monoxide and carbon dioxide
Nitrogen oxides (NO_x)
- **5.3 Advice for firefighters**
- **Protective equipment:**
Wear fully protective suit.

(Contd. on page 3)

GB

Safety Data Sheet

according to UK REACH (SI 2020/1577) as amended

Date of issue: 11.09.2026

Version number: 1.1 (replaces version 1.0)

Revision date: 07.11.2012

Trade name: Fmoc-Gln(trt)-OH >98.5%

(Contd. of page 2)

Wear self-contained respiratory protective device.

Do not inhale explosion gases or combustion gases.

· Additional information

Cool endangered receptacles with water spray.

Collect contaminated fire fighting water separately. It must not enter the sewage system.

*** SECTION 6: Accidental release measures****· 6.1 Personal precautions, protective equipment and emergency procedures**

Avoid formation of dust.

Ensure adequate ventilation

Wear protective equipment. Keep unprotected persons away.

· 6.2 Environmental precautions: No special measures required.**· 6.3 Methods and material for containment and cleaning up:**

Pick up mechanically.

Send for recovery or disposal in suitable receptacles.

· 6.4 Reference to other sections

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

*** SECTION 7: Handling and storage****· 7.1 Precautions for safe handling**

Prevent formation of dust.

Ensure good ventilation/exhaustion at the workplace.

Avoid direct contact (skin/eye contact, ingestion and/or inhalation of fume/mist/dust) with the product in the undiluted form.

· Information about fire - and explosion protection:

There is a risk of a dust explosion if the following conditions are met:

- The substance is present in very finely distributed form (powder, dust).

- The substance is whirled up in sufficient quantity in the air.

- An ignition source is present (flame, spark, electrostatic discharge, etc.)

· 7.2 Conditions for safe storage, including any incompatibilities**· Storage:****· Requirements to be met by storerooms and receptacles:** No special requirements.**· Information about storage in one common storage facility:** Store away from oxidising agents.**· Further information about storage conditions:** Store in cool, dry conditions in well sealed receptacles.**· 7.3 Specific end use(s)** No further relevant information available.*** SECTION 8: Exposure controls/personal protection****· 8.1 Control parameters****· Ingredients with limit values that require monitoring at the workplace:** Not required.**· Additional information:** The lists valid during the making were used as basis.**· 8.2 Exposure controls****· Appropriate engineering controls**

Provide general and/or local exhaust ventilation to control airborne levels below the exposure limits.

If exposure limits have not been established, maintain airborne levels to an acceptable level.

(Contd. on page 4)

— GB —

Safety Data Sheet

according to UK REACH (SI 2020/1577) as amended

Date of issue: 11.09.2026

Version number: 1.1 (replaces version 1.0)

Revision date: 07.11.2012

Trade name: Fmoc-Gln(trt)-OH >98.5%

(Contd. of page 3)

- **Individual protection measures, such as personal protective equipment**
- **General protective and hygienic measures:**
The usual precautionary measures are to be adhered to when handling chemicals.
Do not breathe dust
Avoid contact with the eyes and skin.
Do not eat, drink, smoke or sniff while working.
Wash hands before breaks and at the end of work.
Depending on the degree of exposure, periodic medical examination is suggested.
Do not carry product impregnated cleaning cloths in trouser pockets.
- **Respiratory protection:** Not necessary if room is well-ventilated.
- **Hand protection**



Protective gloves.

Use gloves tested and approved under appropriate government standards such as NIOSH (US) or EN374 (EU).

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

- **Material of gloves**
The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer.
- **Penetration time of glove material**
The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.
- **Eye/face protection**



Safety glasses with side-shields conforming to EN166.

Use equipment for eye protection tested and approved under appropriate government standards such as NIOSH (US) or EN 166(EU).

- **Body protection:**



Protective work clothing

Body protection must be chosen depending on product properties, activity and possible exposure.

* SECTION 9: Physical and chemical properties

- **9.1 Information on basic physical and chemical properties**
- **General Information**
- **Physical state** Solid
- **Colour:** White
- **Odour:** Mild
- **Odour threshold:** Not determined
- **Melting point/freezing point:** 165 °C
- **Boiling point or initial boiling point and boiling range** Undetermined
- **Flammability** Product is not flammable.
- **Lower and upper explosion limit**
- **Lower:** Not determined
- **Upper:** Not determined
- **Flash point:** Not applicable
- **Decomposition temperature:** Not determined
- **pH (0.09 g/l) at 20 °C** 5.4

(Contd. on page 5)

Safety Data Sheet

according to UK REACH (SI 2020/1577) as amended

Date of issue: 11.09.2026

Version number: 1.1 (replaces version 1.0)

Revision date: 07.11.2012

Trade name: Fmoc-Gln(trt)-OH >98.5%

(Contd. of page 4)

· Viscosity:	
· Kinematic viscosity	Not applicable
· Dynamic:	Not applicable
· Solubility	
· water at 20 °C:	0.09 g/l
· Partition coefficient n-octanol/water (log value)	4.2 log POW
· Vapour pressure:	Not applicable
· Density and/or relative density	
· Density:	Not determined
· Relative density	Not determined
· Vapour density	Not applicable
· Particle characteristics	Not determined
· 9.2 Other information	
· Appearance:	
· Form:	Powder
· Important information on protection of health and environment, and on safety.	
· Ignition temperature:	Not determined
· Explosive properties:	Product is not explosive. However, formation of explosive air/dust mixtures are possible.
· Change in condition	
· Evaporation rate	Not applicable.
· Information with regard to physical hazard classes	Not applicable

* SECTION 10: Stability and reactivity

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability**
- **Thermal decomposition / conditions to be avoided:**
There is a risk of a dust explosion if the following conditions are met:
 - The substance is given in very finely distributed form (powder, dust).
 - The substance is whirled up in sufficient quantity in the air.
 - An ignition source is present (flame, spark, electrostatic discharge, etc.)
- **10.3 Possibility of hazardous reactions** Reacts with strong oxidising agents.
- **10.4 Conditions to avoid** Heat and static discharge.
- **10.5 Incompatible materials:** Strong acids and oxidising agents
- **10.6 Hazardous decomposition products:**
Carbon monoxide and carbon dioxide
Nitrogen oxides (NO_x)

* SECTION 11: Toxicological information

- **11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008**
- **Acute toxicity** Based on available data, the classification criteria are not met.
- **Primary irritant effect:**
- **Skin corrosion/irritation** Based on available data, the classification criteria are not met.
- **Serious eye damage/irritation** Based on available data, the classification criteria are not met.
- **Respiratory or skin sensitisation**
May cause an allergic skin reaction.
- **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
- **Carcinogenicity** Based on available data, the classification criteria are not met.

(Contd. on page 6)
GB

Safety Data Sheet

according to UK REACH (SI 2020/1577) as amended

Date of issue: 11.09.2026

Version number: 1.1 (replaces version 1.0)

Revision date: 07.11.2012

Trade name: Fmoc-Gln(trt)-OH >98.5%

(Contd. of page 5)

- **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.
- **Subacute to chronic toxicity:** Prolonged or repeated skin contact may irritate and cause dermatitis.
- **Additional toxicological information:** May cause sensitisation in susceptible persons.
- **11.2 Information on other hazards**
- **Endocrine disrupting properties** Substance is not listed.

SECTION 12: Ecological information

- **12.1 Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **12.2 Persistence and degradability** biodegradable
- **12.3 Bioaccumulative potential** Product is not expected to bioaccumulate.
- **12.4 Mobility in soil** No further relevant information available.
- **12.5 Results of PBT and vPvB assessment**
- **PBT:** Not applicable
- **vPvB:** Not applicable
- **12.6 Endocrine disrupting properties** The product does not contain substances with endocrine disrupting properties.
- **12.7 Other adverse effects**
- **Additional ecological information:**
- **General notes:**
Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

* SECTION 13: Disposal considerations

- **13.1 Waste treatment methods**
- **Recommendation**
Recommended Hierarchy of Controls:
 - Minimise waste;
 - Reuse if not contaminated;
 - Recycle, if possible; or
 - Safe disposal (if all else fails).
 Small amounts may be diluted with plenty of water and washed away. Dispose of bigger amounts in accordance with Local Authority requirements.
Contact waste processors for recycling information.
Used, degraded or contaminated product may be classified as hazardous waste. Anyone classifying hazardous waste and determining its fate must be qualified in accordance with state and international legislation.

· **UK List of Waste**

HP13	Sensitising
HP14	Ecotoxic

- **Uncleaned packaging:**
- **Recommendation:**
Empty contaminated packagings thoroughly. They may be recycled after thorough and proper cleaning.
Disposal must be made according to official regulations.
Containers, even those that are "empty," may contain residues that can develop flammable and/or hazardous vapours upon heating. Do not cut, drill, grind, weld, or perform similar operations on or near empty containers.

(Contd. on page 7)

Safety Data Sheet

according to UK REACH (SI 2020/1577) as amended

Date of issue: 11.09.2026

Version number: 1.1 (replaces version 1.0)

Revision date: 07.11.2012

Trade name: Fmoc-Gln(trt)-OH >98.5%

(Contd. of page 6)

- **Recommended cleansing agents:** Water, if necessary together with cleansing agents.

SECTION 14: Transport information

- | | |
|---|--|
| · 14.1 UN number or ID number
· ADR, ADN, IMDG, IATA | Not applicable |
| · 14.2 UN proper shipping name
· ADR, ADN, IMDG, IATA | Not applicable |
| · 14.3 Transport hazard class(es)
· ADR, ADN, IMDG, IATA
· Class | Not applicable |
| · 14.4 Packing group
· ADR, IMDG, IATA | Not applicable |
| · 14.5 Environmental hazards:
· Marine pollutant: | No |
| · 14.6 Special precautions for user | Not applicable |
| · 14.7 Maritime transport in bulk according to IMO instruments | Not applicable |
| · Transport/Additional information: | Not dangerous according to the above specifications. |
| · UN "Model Regulation": | Not applicable |

SECTION 15: Regulatory information

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
 - **Poisons Act**
 - **Regulated explosives precursors** Substance is not listed.
 - **Regulated poisons** Substance is not listed.
 - **Reportable explosives precursors** Substance is not listed.
 - **Reportable poisons** Substance is not listed.
 - **Per- and polyfluoroalkyl substances (PFAS)** Substance is not listed.
- **Control Of Major Accident Hazards Regulations 2015 (COMAH)**
- **Named dangerous substances - ANNEX I** Substance is not listed.
- **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

* SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.
 This Safety Data Sheet is in compliance with Regulation (EC) No 1907/2006, Article 31 as amended by Regulation (EU) 2020/878.

- **Department issuing SDS:** Product safety department.
- **Date of previous version:** 07.11.2012
- **Version number of previous version:** 1.0
- **Abbreviations and acronyms:**
 UK- REACH: The UK REACH etc. (Amendment etc.) (EU Exit) Regulations 2019 concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals

(Contd. on page 8)

Safety Data Sheet
according to UK REACH (SI 2020/1577) as amended

Date of issue: 11.09.2026

Version number: 1.1 (replaces version 1.0)

Revision date: 07.11.2012

Trade name: Fmoc-Gln(trt)-OH >98.5%

(Contd. of page 7)

GB- CLP: Retained GB CLP Regulation (EU) No. 1272/2008 as amended for Great Britain on the Classification, Labelling and Packaging of Chemicals.

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

Skin Sens. 1: Skin sensitisation – Category 1

Aquatic Chronic 4: Hazardous to the aquatic environment - long-term aquatic hazard – Category 4

· *** Data compared to the previous version altered.**

GB
